

Technical Information

Primax PT- D

X-ray Film Tank Developer (powder)

Description :

- X - ray developer for (manual) tank processing of all customary X - ray films
- powdery, two parts

Processing Conditions :

developing temperature :	18 °C	developing time :	5 - 6 min.
	20 °C		4 - 5 min.
	22 °C		3 - 4 min.
	24 °C		2 - 3 min.

Replenishment :

For processing of 1 m² X - ray film **300 to 400 ml** developer solution are necessary at standard conditions.

Mix Instructions :

Water (30 - 40 °C)	Part A	Part B	fill up with Water to
 about ¾ of tank	 First dissolve part A completely in water, afterwards add part B under permanent stirring, then fill up to the final volume and mix well	 then fill up to the final volume and mix well	 10,00 l
7,5 l			22,5 l
17,0 l			

Assortment :

Article Number	Product Name	Mix Size for	Capacity for ... m ² film
209133	Primax PT-D	10 x 10 l	240 - 330
209196		5 x 22,5 l	280 – 375

Storage and Handling :

- Until the preparation PT-D has to be kept in the closed original package and stored at temperatures of 10 to 25 ° C.
- All users of photographic processing solutions should exercise the greatest care to avoid the chemicals contacting the skin or eyes . Always wear solution resistant gloves.
- The disposal of used photographic processing solutions has to be occurred in accordance with the legal regulations.