

## Front apron **PRIMAX** soft, Standard style

The Front apron Primax *soft* is normally worn as front protection for short and occasional periods during conventional X-ray.

- appropriate radiation protection for medical personnel
- 0.25 mm Pb, 0.35 mm Pb or 0.50 mm LE
- can be worn in combination with protective neck cover
- convenient Velcro-fastener crosswise closure on hips, therefore easy to open and easy to handle
- sleeveless, with protection of shoulder and part protection of back, with stitched on breast pocket
- comfort of the user is guaranteed by shoulder-padding
- available in different sizes and lengths
- standard colour – Royal blue
- available in five additional colours



Size (L x W)	LE	Article-No.	LE	Article-No.	LE	Article-No.
100 x 60 cm	0,25 mm	224001	0,35 mm	224011	0,50 mm	224031
110 x 60 cm	0,25 mm	224002	0,35 mm	224012	0,50 mm	224032
120 x 60 cm	0,25 mm	224003	0,35 mm	224013	0,50 mm	224033
120 x 75 cm	0,25 mm	224004	0,35 mm	224014	0,50 mm	224034