

## Protection for dental medicine PRIMAX<sup>®</sup> soft

### Dental apron

- 35 mm Pb or 0.50 mm Pb equivalent protection
- for use with patients undergoing dental intra-oral x-ray exposures
- full coverage of front body and shoulder of the patient
- Velcro-fastener in the neck, easy to handle
- available in different sizes and lengths
- standard colour – Royal blue
- all items available in five additional colours



GREEN PETROL DARK BLUE VIOLET RED

	Size	LE	Article-No.
Adult	(50 x 76 cm)	0,35 mm	224500
Adult	(50 x 76 cm)	0,50 mm	224504
Children	(45 x 60 cm)	0,35 mm	224501
Children	(45 x 60 cm)	0,50 mm	224506

### Panoramic dental apron

- with vertebral protection
- 35 mm Pb or 0.50 mm Pb equivalent protection
- for use with patients undergoing panorama dental x-ray exposures and remote x-ray exposures
- individual fixing by three horizontal Velcro closures in the front part, shoulder padding
- available in different sizes and lengths



	Size	LE	Article-No.
Adult	(53 x 67 cm)	0,35 mm	224520
Adult	(53 x 67 cm)	0,50 mm	224524
Children	(46 x 56 cm)	0,35 mm	224521
Children	(45 x 56 cm)	0,50 mm	224525

### Neck cover (with extension)

- used for protection of thyroid and sternum
- for use with patients undergoing dental intra-oral x-ray exposures
- easy to handle neck Velcro-fastener
- neck cover suitable for size 38 – 46 , collar height 6 cms



	Size	LE	Article-No.
	One size (38-46)	0,50 mm	224402