

## Front apron **PRIMAX** soft, Surgeon style

The surgeon apron *Primax soft* is designed to be worn during operating room procedures under the sterile gown. The special wrap around straps can be easily released allowing the apron to fall free without removing or disturbing the sterile gown.

- appropriate radiation protection for medical personnel
- 0.25 mm Pb, 0.35 mm Pb or 0.50 mm LE
- can be worn in combination with protective neck cover
- apron fastening from front body part by two straps going over shoulder, crossing on the back and then going partly to the front again
- Velcro-fastening system and shoulder padding for user's comfort
- breast pocket for personal dosimeter
- available in different sizes and lengths
- available in five additional colours



GREEN PETROL DARK BLUE VIOLET RED



Size (L x W)	LE	Article-No.	LE	Article-No.	LE	Article-No.
100 x 60 cm	0,25 mm	224101	0,35 mm	224111	0,50 mm	224131
110 x 60 cm	0,25 mm	224102	0,35 mm	224112	0,50 mm	224132
120 x 60 cm	0,25 mm	224103	0,35 mm	224113	0,50 mm	224133
120 x 75 cm	0,25 mm	224104	0,35 mm	224114	0,50 mm	224134