

Protective vest and skirt PRIMAX[®] soft

The two-piece wrap around protection, consisting of vest and skirt is an ideal selection for radiological procedures which require personal protection for extended period. This combination guarantees more mobility with division of weight on shoulder and hips.

Protective Vest Primax soft

- weight is optimally divided between the shoulders and hips
- 0.25 or 0.35 mm LE (back) / 0.50 mm LE (front)
- individual fixing by three horizontal Velcro closures in the front part, shoulder padding
- available in different sizes and lengths
- standard colour – Royal blue
- available in five additional colours



VEST

Size (L x W)	LE (front/back)	Article-No.	LE (front/back)	Article-No.	LE (front/back)	Article-No.
S (60 x 85-94 cm)	0,50/0,25 mm	224315	0,50/0,35 mm	224324	0,50/0,50 mm	224329
M (60 x 95-104 cm)	0,50/0,25 mm	224316	0,50/0,35 mm	224318	0,50/0,50 mm	224327
L (60 x 105-114 cm)	0,50/0,25 mm	224317	0,50/0,35 mm	224323	0,50/0,50 mm	224326
XL(65 x 115-124 cm)	0,50/0,25 mm	224319	0,50/0,35 mm	224325	0,50/0,50 mm	224328

Protective vest and skirt PRIMAX[®] soft



Protective Skirt Primax soft

- 0.25 mm; 0.35 mm or 0.50 mm LE (back) / 0.50 mm LE (front)
- easy skirt handling by wide double waist closing, Velcro closure and belt with quick release buckle
- available in different sizes and lengths
- standard colour – Royal blue
- available in five additional colours



GREEN PETROL DARK BLUE VIOLET RED

SKIRT

Size (L x W)	LE (front/back)	Article-No.	LE (front/back)	Article-No.	LE (front/back)	Article-No.
S (60 x 85-94 cm)	0,50/0,25 mm	224335	0,50/0,35 mm	224341	0,50/0,50 mm	224344
M (60 x 95-104 cm)	0,50/0,25 mm	224336	0,50/0,35 mm	224339	0,50/0,50 mm	224345
L (60 x 105-114 cm)	0,50/0,25 mm	224337	0,50/0,35 mm	224340	0,50/0,50 mm	224346
XL(65 x 115-124 cm)	0,50/0,25 mm	224338	0,50/0,35 mm	224342	0,50/0,50 mm	224347