

# Wofasept® FL spray disinfection

Ready-to-use surface disinfectant



For odourless disinfection of washable surfaces in the working area of medical practices, hospitals and other medical facilities as well as laboratories, food industry and gastronomy. Can be used when alcoholic disinfection is undesirable.

**New:**  
Bottle with  
practical foam  
and spray  
function



## Highly effective microbicide

Wofasept® FL spray disinfection effectively kills bacteria, funguses, lipophilic viruses (including HBV/HIV) and human rotaviruses



## No protein related faults

Insensitive to contamination such as blood, protein, stool, sputum, pus, etc.



## No allergic potential

Wofasept® FL spray disinfection is free of chlorine, aldehyde and phenol and does not contain any substances from which an allergic potential is to be expected.



## Application and dosage:

Spray or foam undiluted wipe and wipe away with a dry tissue after application, if required.

	Application time
for surface disinfection in hospitals, practices and other medical facilities	
• kills bacteria and funguses	5 min
• kills hepatitis B and HIV	1 min
• kills herpes viruses	15 min
• kills human rotaviruses	1 min
• kills vacciniaviruses	5 min
• kills noroviruses	30 min
• kills influenza A-viruses	1 min
Disinfects solariums and washable surfaces	5 min
for food industry	30 min



