

Double sided apron **PRIMAX** soft

The coat apron Primax *soft* provides an appropriate front and back protection in full length. Preferable it will be used during X-ray operations in cardiology, angiography and endoscopy and where required also for assisting personnel in the operating theatre.

- perfect as all-around protection for all examinations
- 0.25 mm LE, 0.35 mm LE or 0.50 mm LE (front) / 0.25 mm LE (back)
- fixed belt and shoulder padding to reduce the pressure on the shoulders
- care should be taken to ensure adequate overlapping when fastening (min. 10 cm)
- can be worn in combination with protective neck cover
- available in different sizes and lengths
- standard colour – Royal blue, other colours on request



Size (L x W)	LE (front/back)	Article-No.	LE (front/back)	Article-No.	LE (front/back)	Article-No.
100 x 60 cm	0,25/0,25 mm	224201	0,35/0,25 mm	224211	0,50/0,25 mm	224221
110 x 60 cm	0,25/0,25 mm	224202	0,35/0,25 mm	224212	0,50/0,25 mm	224222
120 x 60 cm	0,25/0,25 mm	224203	0,35/0,25 mm	224213	0,50/0,25 mm	224223
120 x 75 cm	0,25/0,25 mm	224204	0,35/0,25 mm	224214	0,50/0,25 mm	224224