

Gonad shielding PRIMAX[®] soft

Half apron

- for protection of patients gonads during radioscopy of thorax and upper abdomen
- 0.50 mm Pb or 1.00 mm LE
- waist fastening by simple buckle
- available in 4 sizes (for toddlers, children, adults)
- standard colour - Royal blue
- available in five additional colours



Size	LE	Article-No.
Children 1 (20 x 25 cm)	0,50 mm	224410
Children 2 (30 x 30 cm)	0,50 mm	224411
Adult 1 (37 x 40 cm)	0,50 mm	224412
Adult 2 (45 x 50 cm)	0,50 mm	224413
Children 1 (20 x 25 cm)	1,00 mm	224415
Children 2 (30 x 30 cm)	1,00 mm	224416
Adult 1 (37 x 40 cm)	1,00 mm	224417
Adult 2 (45 x 50 cm)	1,00 mm	224418

Ovary Shield

- 1.00 mm LE, set of 2 parts
- protection of the ovaries of female patients during X-ray examinations
- available in 2 sizes

Medium (9 x 12 cm)	1,00 mm LE	224430
Large (11 x 18 cm)	1,00 mm LE	224431
Set of both sizes	1,00 mm LE	224432



Scrotum shields

- 1.00 mm LE, available in 2 sizes
- protection of the ovaries of male patients during X-ray examinations

Medium (9 x 12 cm)	1,00 mm LE	224435
Large (11 x 18 cm)	1,00 mm LE	224436
Set of both sizes	1,00 mm LE	224437



Flexible gonad protection (Shadow shields)

- 1.00 mm LE , set of 3 pieces
- used when scrotum and ovary shields cannot be applied

Set (10x12/14x17/18x21 cm)	1,00 mm LE	224440
----------------------------	------------	--------

